

Logic and Accuracy Testing Checklist



ClearAccess and ClearCast

Precinct / Vote Center	
-------------------------------	--

ClearAccess

ITEM	DESCRIPTION	COMPLETED
1	Load Election using ADF on USB Drive	
2	Uncheck print stubs in settings	
3	Test Vote -Using all available inputs	
4	Duplicate Ballot(s) from the test deck	
5	IMPORTANT: Set Options to print ballots without stubs	
6	Set for Election	
7	Shut Down	
8	Remove USB Drive	



ClearCast

ITEM	DESCRIPTION	COMPLETED
1	Insert 2 USB sticks into Cast Unit and plug it in	
2	Verify ClearCast Name from Power on Report match precinct	
3	Login as Maintenance	
4	IMPORTANT: Calibrate Scanner	
5	Login as System Administrator	
6	Load Election and Set Correct Precinct	
7	Set for Logic and Accuracy	
8	Scan Test Deck including duplicated ballot(s) from ClearAccess	
9	Stop scanning (touch ClearCast icon in upper left)	
10	Display and print Summary Results Tape	
11	Leave Test Ballots on Cast	
12	Compare results to expected results; collect ballots and tape.	
13	End Test and Shut Down	
14	Remove USB Drives for Central Merge	
	After Central Merge and w/cleaned USB Drives	
15	Login as System Administrator	
16	Load Election and Set Correct Precinct	
17	Set for Election	
18	Shut down	